

**Piezo Mini**

# MINI DENTAL PIEZO SURGERY

Minisize and lightweight  
dual-screen design for smart dual control



# MINI DESIGN

## COMPACT & EFFICIENT

Balancing minimally invasive techniques & power output

# 5"

MAIN SCREEN



Full-touch interface for intuitive operation

# 1.6"

ROUND SCREEN



Real-time data display and adjust

### OPTIMIZED DIMENSIONS

perfectly designed for a streamlined workspace

### SLIM HANDPIECE

Operates at full power for 10 mins with a handpiece temperature of only 33°C, ensuring safety and comfort



Bone Mode



Endo Mode



Perio Mode



Custom Mode



Clean Mode

### PARAMETER

<b>VOLTAGE</b>	100~240V 50/60Hz 120VA	<b>RIDAL AMPLITUDE</b>	20-200µm
<b>FUSES</b>	2×F10A 250V	<b>LATERAL AMPLITUDE</b>	≤30µm
<b>FREQUENCY</b>	26~38Hz	<b>WATER FLOW</b>	60-100ml